

# AGRICULTURAL INDUSTRY

## Soils

Filtration and microfiltration products for applications in agricultural analysis.

The following filter paper and microfiltration products are ideal for soil analysis applications.



PROCESS	TECHNIQUE	TYPE OF FILTER	REF. CHM
Phosphorus determination	Filtration (Funnel/Buchner Funnel)	Quantitative filter paper (Nitrogen less and free of phosphates)	Under request
Nitrogen determination		Weighing boats	E4002
Soluble sulphates determination		Quantitative filter paper	F2044-
Carbonates determination		Quantitative filter paper	F2040-

